



PORTABLE/FOLDING POWER WHEELCHAIR COMPARISON CHART

Model	Weight (Without Battery)	Top Speed	Range per Charge	Turning Radius	Weight Capacity	Folding Type	Battery Type	Airline-Friendly
Aerolux	29 lbs	3.7 mph	11.2 miles	26"	300 lbs	Cross-Folding	Lithium-Ion	Yes
Dash Ultra-Lite	39 lbs	4.3 mph	9.5 miles	30"	265 lbs	Manual Fold	Lithium-Ion	Yes
Feather Power Chair	33 lbs	4 mph	13 miles	35"	350 lbs	Manual Fold	Lithium-Ion	Yes
Go Chair	35 lbs	3.8 mph	11.2 miles	25.5"	300 lbs	Disassembles	Lithium-Ion	Yes
Golden Ally	37 lbs	4 mph	9 miles	28.5"	330 lbs	Manual Fold	Lithium-Ion	Yes
Golden Ally Pro	37 lbs	5 mph	13 miles	24"	330 lbs	Manual Fold	Lithium-Ion	Yes
Golden Ally Sport	35 lbs	4 mph	9.5 miles	28"	265 lbs	Manual Fold	Lithium-Ion	Yes
Golden BuzzAbout	34.8 lbs (heaviest piece)	4 mph	11 miles	28.5"	300 lbs	Disassembles	Lithium-Ion	Yes
Buzzaround CarryOn	46 lbs (heaviest piece)	3.85 mph	12.5 miles	30"	300 lbs	Manual Fold	Lithium-Ion	Yes
Jazzy Carbon	43.6 lbs	3.7 mph	9.3 miles	24"	300 lbs	Manual Fold	Lithium-Ion	Yes
Jazzy Ultra Light	33 lbs	3.7 mph	16.4 miles	24"	300 lbs	Manual Fold	Lithium-Ion	Yes
LiteRider Envy	37 lbs	3.5 mph	15.5 miles	28.5"	300 lbs	Disassembles	SLA	No
Travel Buggy City 2 Plus	40 lbs (heaviest piece)	4.5 mph	15 miles	32"	395 lbs	Quick Fold	Lithium-Ion	Yes
WHILL Model C2	42 lbs (heaviest piece)	5 mph	11 miles	29.9"	300 lbs	Disassembles	Lithium-Ion	Yes
WHILL Model F	52.9 lbs	3.7 mph	12.4 miles	30.7"	250 lbs	Manual Fold	Lithium-Ion	Yes
Zoomer Power Chair	50 lbs	3.7 mph	8 miles	29"	275 lbs	Manual Fold	Lithium-Ion	Yes

Last Updated: May 2025 - Please note that models shown are Freedom Mobility Solutions Models ONLY.