



# FEATURED HI-LOW BEDS



## Full Electric Standard/Low Bed

- Bariatric size available (not compatible with Halo).
- Adjustable height with dual casters.
- Easy delivery with split spring design.
- Independent operation of head, foot, and hi/lo motors.
- **Also available for rent!**

Halo Compatible



## Prime Care Bed Model P903

- Easy expansion to 42" bariatric width.
- Ultra-low 7" deck height for fall risk prevention.
- Dual pedal floor lock system keeps bed securely in place.
- The Backlit Hand Pendant includes Auto Contour, Cardiac Comfort Chair, and under-bed night light.

Halo Compatible



## SonderCare - Aura Hospital Bed

- Ultra-low 10" bed height.
- Infinite positioning.
- High height of 32" allows caregivers to work at a comfortable height.
- Available in 39" and 48" widths.
- Comes with a washable incontinence cover.

Halo Compatible



9"H Optional Fall Mat stows under bed.

## Prius Descend Ultra-Low Bed

- Enhanced safety with ultra-low floor design for fall prevention.
- Electric height adjustment from 2.7" to 25".
- Customizable backrest, footrest, and Trendelenburg functions.
- Battery backup for emergency lowering.
- Auto-lock for added safety.



## Supernal 3 Bed

- Available in Twin, Full, Queen, or Dual King sizes.
- Adjustable bed deck range of 12"-24".
- Aesthetically pleasing homestyle design.
- Compatible with a patient lift (with added casters).
- Shown with Wonder Pole, a Halo alternative.



## Advantage ReadyWide Bed

- ReadyWide two-click expansion to 39" or 42" width.
- Comfort Chair, Trendelenburg, and Reverse Trendelenburg positions.
- Hi-Lo adjustment from 7-7/8" to 30".
- For an upgrade, check out the Encore ReadyWide Bed.

BROWSE ALL



HOSPITAL BEDS  
& MATTRESSES

- ✓ **Trained Bed and Mattress Experts to Assist You**
- ✓ **2 Showrooms to Test and Try Options**
- ✓ **White Glove Delivery and Setup**
- ✓ **Freedom Mobility Promise with Every Purchase**



📍 Hanover, MD Showroom



📞 +1 (888) 414-6553

📍 Now Serving MD, DC, and NOVA

🌐 [FreedomMobilitySolutions.com](http://FreedomMobilitySolutions.com)



# FEATURED MATTRESSES



## SPAN Geo-Mattress UltraMax

- Advanced foam mattress with three patented features, including Shear Transfer Zones for preventing pressure points
- Bi-directional stretch cover, fluid-proof top cover, and over 800 articulating cells for enhanced comfort and durability



## Prius VLX Foam Mattress

- Tri-layered foam mattress with cool visco elastic top layer
- Zoned comfort foam, heel soft foam, and firm support base
- Full-length firm side rails, urethane 4-way stretch top cover
- Available in standard, bariatric, and custom sizes.



## Prius Salute RDX Mattress System

- Budget-friendly premium quality alternating air mattress
- Therapeutic pressure
- Power outage inflation capability
- 3-in-1 alternating air cycle, 8 adjustable comfort levels
- Low friction, low shear cover for maximum protection and patient comfort
- **Also available for RENT**



## Span America CustomCare Convertible Mattress

- Lateral rotation and alternating pressure in pump operation
- Foam layer on top for enhanced comfort
- Self-enclosed vertical air cylinders for pump-free use
- Easily converts to a powered surface with a convenient control unit
- Shear transfer zone technology for skin protection

Top Name Brands



**Make an Appointment to Come Visit  
One of Our Two Showrooms**

### MARYLAND

7271 Park Circle Drive  
Unit 100  
Hanover, MD 21076

### VIRGINIA

21055 Ashburn Crossing Dr.  
Unit 120  
Ashburn, VA 20147