

# Supernal™

SLEEP SYSTEM

The **Supernal Sleep System** is a one-stop shop for all of your mattress needs. The two mattresses in this line are the perfect matches for our **Supernal Recliner Plus** and **Supernal Hi-Low beds**, as well as an outstanding choice for replacing an existing adjustable bed.

Two mattress styles to meet a wide range of needs:

## Ascent I and the Soft Touch

### Ascent I Elements:

-  Includes state-of-the-art progressive laminated foams to provide a healthy skin environment.
-  Border edge protection.
-  Designed with sleep technologies that promote excellent comfort and healthy bedding solutions.
-  Manufactured in an FDA approved facility that has decades of manufacturing experience to ensure the highest quality.

### Soft Touch Elements:

-  Core made of cool, form-fitting memory foam that conforms to the shape of the person who is using it.

-  Border edge protection.
-  Memory foam users credit their wholesome night's rest to the absolute comfort of the soft memory foam.

Two Mattress Cover Options:  
**Natural Bamboo Quilted Cover** and the **Stretch Fabric Cover**

### Natural Bamboo Quilted Cover Elements:

-  Comfortable and fresh sleeping environment.
-  Soft, luster, and silky feel.
-  Natural properties: anti-bacterial, deodorizing, and breathable.
-  Comprised of a Joma wool filling, which absorbs and causes moisture to evaporate, for maximum comfort.

### Stretch Fabric Cover Elements:

-  A unique, waterproof cover that feels soft, just like fabric.
-  Can be cleaned with household cleansers.
-  Will protect your mattress for many years.



SUPERNAL RECLINER PLUS

# INTRODUCING: Supernal™ RECLINER BEDS



Currently seeking distributors worldwide.

All images are the property of Transfer Master Products, Inc.

505 W. Williams PO Box 917 Postville, IA 52162

Call: 563-864-7364 or 800-475-8122

TransferMaster.com • E-mail: Sales@TransferMaster.com



# Introducing Supernal<sup>™</sup> RECLINER BEDS

The **Supernal Recliner Bed System** is designed specifically for individuals who need hospital bed features at home but do not want the hospital room look. The bed is upholstered for maximum comfort and designed in a way that conceals all electrical components, giving it the appearance of an ordinary bed. The **Supernal Recliner Plus** bed enables its users to proudly receive visitors in their bedrooms without sacrificing their dignity or compromising on functionality.

The patent-pending bed frame combines the quality, elegance, craftsmanship, and cost-effectiveness of American, European, and Asian engineering. The proprietary Supernal Sleep mattresses are designed and manufactured exclusively in the United States.

The **Supernal Recliner Plus** is both high-quality and affordable. Its features include:

- ☞ European-style head tilt to cradle pillow
- ☞ Wall hugging; night table is always within reach of the user
- ☞ Head and foot adjustability

- ☞ Illuminated wireless control; DC-powered with emergency power-down
- ☞ Compatible with 110 or 220 volts, robust massage motors
- ☞ Low settings for easy transfer into and out of the bed
- ☞ Available in multiple sizes, including Twin, Full, Queen, and Dual King



The **Supernal Hi-Low** has all of the features of the Recliner Plus model, with the following specific capabilities:

- ☞ 10" height adjustment
- ☞ Hi-Low "bolt on" for easy delivery, set-up, and storage
- ☞ Low settings for easy transfer into and out of the bed

See the back cover of the brochure for detailed information about the Supernal Sleep System mattresses.

Supernal<sup>™</sup>  
RECLINER PLUS

HEAD + FOOT  
ADJUSTABLE

Supernal<sup>™</sup>  
HI-LOW

HEAD + FOOT  
HI-LOW ADJUSTABLE



Bed Flat

YES

YES



← Head Tilt

Head Tilt

YES

YES



Head Fully Raised

YES

YES



Head and Foot Raised

YES

YES



Hi-Low Option

10 Inches

NO

YES