


THE HALO SAFETY RING



PRODUCT FEATURES

The Halo Safety Ring is the only patented institutional bed mobility device that prevents lateral mattress movement and gap space in the most critical zones: 1, 2, 3, & 4. The Halo Safety Ring meets and exceeds the HBSW and FDA's final dimensional guidelines against bed entrapment. The Halo Safety Ring passes the new "Bed System Measurement Device" pass/fail test that has been approved for Bed System testing by the FDA and state surveying teams.


- + Provides exit points for neck, head, & chest and eliminates the four most critical zones of entrapment.
- + Meets FDA guidelines for patient transfers & therapy bed side assistant for both patient transfers & therapy.
- + Promotes bed mobility & increases environmental independence.



- **Zone 1:** Within the rail.
- **Zone 2:** Between the top of the compressed mattress and the bottom of the rail, between rails ports.
- **Zone 3:** Between the rail and the mattress.
- **Zone 4:** Between the top of the compressed mattress and the bottom of the rail, at the end of the rail.
- **Zone 5:** Between the split bed rails.
- **Zone 6:** Between the end of the rail and the side edge of the head or foot board.
- **Zone 7:** Between the head or footboard and the mattress end.

BED FRAME MOUNT

MOUNTS TO BED FRAME: Both Halo ALC and Institutional Bed Model mount to bed frame for stability.



BED HANDLE ROTATES

ROTATING HANDLE:

Rotates 90° either direction to accommodate different user positions and allow ease of bed exit.



HALO SAFETY RING INSTITUTIONAL BED MODEL



FEATURES

- Available in single and double sided for institutional style beds.
- Easy to install turn-buckle mounting bracket.
- 1000 lb. weight capacity steel construction

HALO SAFETY RING INSTITUTIONAL BED MODEL

MODEL NUMBER	DESCRIPTION
77101	One sided circle device with mattress stay bar and hardware.
77102	Both sides circle device with hardware.
77103	One sided circle device with hardware for carrol.
77113	Halo Safety Ring, adjustable length bed (ALB) unit complete.

HALO ALC ASSISTED LIVING/HOME CARE MODEL



FEATURES

- The universal width and height adjustment accommodates all bed styles (up to queen) and allows bed to be off the wall, preventing bed roll-outs & falls.
- Comes standard with a universal height adjustment system to accommodate all bed systems' mattress and box spring heights.
- 500 lb. weight capacity steel construction.

HALO SAFETY RING ASSISTED LIVING / HOMECARE MODEL

MODEL NUMBER	DESCRIPTION
77121	Halo ALC (Assisted Living Centers) complete kit.

GAP SPACE ELIMINATOR




FEATURES

- Eliminates Zone 7 gap space.
- Designed to eliminate gaps between bed and mattress in increments of 2", 4", and 6".
- Units may be used individually or together, depending on size of gap to be eliminated.
- Each unit has a removable machine washable cover and a three strap system.

GAP SPACE ELIMINATOR

MODEL NUMBER
77117

CIRCLE PAD PROTECTOR



FEATURES

- Prevents bruising from any medical disorder that might cause the skin to become thin, frail, or to bruise.
- Promotes increased safety and a better sleeping environment.
- Installs with simple hook and loop attachments.

CIRCLE PAD PROTECTOR

MODEL NUMBER
77118