



BED RAIL OPTION

Available with Bed Package (PB42BARBED-PKG) or sold separately (#PB8035).



THERAPEUTIC SUPPORT SURFACES

Bariatric Foam Mattress (#PBHB42MAT).

Bariatric Bed

HEAVY-DUTY, FULL-ELECTRIC

The **ProBasics Bariatric Homecare Bed** is a heavy-duty, full-electric bed frame capable of supporting patients up to 600 lbs. (total distributed weight capacity = 700 lbs.). Designed with the bariatric patient in mind, the bed frame extends to the edges of the bed deck for better support than most other split-spring bariatric beds. The hand pendant features dynamic client positioning, allowing the consumer or caregiver to change the positioning of the head and foot sections as well as the bed height. Despite its size, bed setup within the home is easy - no special tools are needed.

FEATURES

- Heavy-duty frame ensures strength and patient safety.
- Larger sleep surface than conventional bed.
- Split-pan design with removable bed ends is easy to setup.
- Motorized positioning of upper body and knees as well as height adjustment provides maximum comfort and versatility to patient and caregiver.
- High impact bed end panels are durable, impact and scratch resistant and washable for easy cleaning.
- Hand control allows for multiple bed positioning for patients.
- Two locking and two non-locking casters for secure positioning.
- Emergency crank provided in the event of a power failure.
- HCPCS Code: E0301*

PRODUCT SPECIFICATIONS

| | |
|---------------------------------|---|
| Bed Height (Floor to Bed Deck): | 18.5" min. – 26" max. |
| Sleep Surface: | 80" L x 42" W |
| Weight Capacity: | Total weight capacity 700 lbs Max patient weight 600 lbs |
| Assembled Bed Weight: | 276 lb |
| Limited Warranty: | 5 year on frame 1 year on electrical components |

Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce the risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment, contact your Compass Health Brands' sales representative for a copy of our "Bed Rail Safety Guide" or visit US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

| Item # | Description | UOM |
|----------------|---|------|
| PB42BARBED | ProBasics Bariatric 42" Full Electric Bed | each |
| PB42BARBED-PKG | ProBasics Bariatric 42" Full Electric with Foam Mattress and Half-Length Brown Vein Rails | each |

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.