

# SEATA

Ergonomic Rollator



2023 Winner!  
The Chicago Athenaeum  
Museum of Architecture and Design



*Sit in Comfort : Walk with Confidence*

# ORDER NOW!

Late March Delivery



**STRONGBACK**  
mobility 

# SEATA

Ergonomic Rollator

## Features



### UNIQUE

Full-size Seat and Ergonomic Backrest for Unmatched Comfort



### INNOVATIVE

Our Design Enables a Stable and Upright Walking Position




### AWARD-WINNING

Lightweight and Foldable Design for Easy Handling

#### Measurements and Weights

Model	SEATA Rollator	
Seat Width	19 in	48 cm
Seat Depth	15 in	38 cm
Seat Height	20 in	50 cm
Seat Angle	12°	
Front Wheel Size	9 in	23 cm
Rear Wheel Size	12 in	29 cm
Overall Width	27 in	68 cm
Overall Length	30 in	77 cm
Overall Height	30-39 in	79-99 cm
Backrest Height	19 in	48 cm
Backrest Angle	25°	
Size Folded	30 x 10 x 31 in	77 x 25 x 80 cm
Total Weight	20 lb	9.3 kg
Maximum Load	260 lb	118 kg

## Sit in Comfort : Walk with Confidence



**Strongback Mobility USA, LLC**

Virginia Beach, VA 23454 • (757) 349-7790

usa@strongbackmobility.com • www.strongbackmobility.com