

# BESTSTAND™

## SA500 | SIT-TO-STAND LIFT

The BestStand™ SA500 is one of the latest addition to our ProCare Sit-to-Stand Solutions. Its ergonomic design ensures safe and smooth patient sit-to-stand transfer.




- **500lb** Capacity
- Full electric with optional spare battery
- Wide Lifting range
- Adjustable knee pad
- Compact base for easy maneuver
- Detachable foot plate

# | DIMENSIONS


Key dimensions of SA500	
maximum capacity	500Lb / 228 Kg
external base width (closed)	26.9"
internal leg width (opened)	35.3"
foot plate clearance	13.8"
upright mast height	41.7"
overall base length	41.5"
base height	5"
minimum lifting height with spreader bar	36.1"
maximum lifting height with spreader bar	66.5"
knee pad length extension range	2.8"
maximum boom width	18.2"
knee pad height extension range	20.2-23.7"
net weight	116lb / 52.6kgs
shipping weight	131lb / 59.4kgs

\* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

## | SA500



TOP VIEW OF SA500



SIDE VIEW OF SA500