

BESTLIFT™

PL500 | FULL BODY PATIENT LIFT

The BestLift™ PL500 is one of the latest ProCare products. Ergonomically designed to reduce and eliminate back injuries among professional caregivers. Providing comfort and dignity for users.




- **500lb** capacity
- Lift from floor and over high beds
- 6-point spreader bar for extra comfort
- Full electric lift with optional spare battery
- Optional scale.

| DIMENSIONS


Key dimensions of PL500	
maximum capacity	500Lb / 228 Kg
external base width (closed)	26.9"
internal leg width (opened)	38.5"
upright mast height	50.6"
overall base length	47.6"
base height	4.9"
minimum lifting height with spreader bar	28.7"
maximum lifting height with spreader bar	72.2"
net weight	98lb / 44.5kgs
shipping weight	122lb / 55.3kgs

* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

| PL500



TOP VIEW OF PL500



SIDE VIEW OF PL500