



# PREMIUM ROLLATOR COMPARISON CHART

Model	Size	Weight	Frame	Wheels	Seat Width	Seat Height	Handle Height	Weight Capacity	Best For
Carbon Ultralight	Compact	10.6 lbs	Carbon Fiber	Puncture-proof	16"	22.4"	30.7" to 35.4"	285 lbs	Travel and tight spaces
Trillium Rollator	Petite / Mini	14 lbs	Aluminum	Polyurethane	18"	18"	29" to 32"	300 lbs	Comfort and everyday use
Xpresso Zero	Petite / Mini	16 lbs	Aluminum	High-traction rubber	17.5"	18"	28.5" to 31.5"	300 lbs	Stability and control
Carbon Ultralight	Regular	10.6 lbs	Carbon Fiber	Puncture-proof	16.5"	24"	31.5" to 37.4"	285 lbs	Balanced everyday mobility
Trillium Rollator	Regular	14 lbs	Aluminum	Polyurethane	18"	21"	31.5" to 34.5"	300 lbs	Comfort and support
Xpresso Zero	Regular	16 lbs	Aluminum	High-traction rubber	17.5"	21"	31" to 35"	300 lbs	Indoor and outdoor use
Carbon Ultralight	Wide	10.6 lbs	Carbon Fiber	Puncture-proof	18.5"	24"	31.5" to 37.4"	285 lbs	Extra room and comfort
Trillium Rollator	Tall	14 lbs	Aluminum	Polyurethane	18"	24"	34.5" to 37.5"	350 lbs	Taller users
Xpresso Zero	Tall	17 lbs	Aluminum	High-traction rubber	17.5"	24"	33" to 37"	350 lbs	Maximum stability