

2-Zone & 3-Position Lift Chair Comparison Chart

Freedom Mobility Solutions | Simple and Versatile Lift Chair Positioning Options

Model	Model #	Positioning	Sizes	Weight Capacity	Heat / Massage	Power Headrest / Lumbar	Charging	Other Features
Golden Cloud Classic	PR510	2-Zone MaxiComfort	S/M, M/L, M/W	400-600 lbs.	Optional	No / No	M/W only	MaxiComfort positioning; custom bucket seat; plush biscuit back
Golden MaxiComforter	PR535	2-Zone MaxiComfort	6 sizes	400-600 lbs.	Optional	No / No	No	Petite through wide sizing; signature seat; optional footrest extension
Golden Comforter	PR531	3-Position / 1 Zone	7 sizes	400-500 lbs.	Optional	No / No	USB - hand control	Signature seat; pocketed coil seat; optional footrest extension
Pride Genesis	LC-150	3-Position	S / M	325 lbs.	No	No / No	No	Fully padded chaise; magnetic hand-control holders; battery backup
Pride Evolution	LC-435	Enhanced 3-Position*	PW / S / M / L / XL / XXL	400-600 lbs.	No	No / No	USB-A/C or USB-C*	Cup holder and storage on select sizes; fully padded chaise

Note: Golden's Comfort Zone terminology and traditional 3-position terminology describe different positioning systems. A 2-Zone chair provides two independently adjustable comfort areas, while a 3-position chair moves through upright, partial recline, and full recline positions. *Pride Evolution positioning and charging features vary by size.

FOR FREEDOM MOBILITY SOLUTIONS EXCLUSIVE USE ONLY