

Petite Lift Chair Comparison Chart

Freedom Mobility Solutions | Petite and Petite-Wide Lift Chair Models

Model	Petite Size	Recommended Height	Weight Capacity	Seat Width	Seat Depth	Seat-to-Floor	Positioning	Heat / Massage	Power Headrest / Lumbar	Charging	Other Features	Wall Clearance
Golden Comforter PR531	Petite/Small	Not specified	375 lbs.	21"	19"	19"	3-Position 1 Comfort Zone	Optional Heat & Vibrating Massage	No / No	USB - hand control	Signature seat; pocketed coil seat; SmartTek diagnostic system; stain defense; optional footrest extension; 9V battery backup	16"
Pride Evolution LC-435	Petite Wide	5'3" & under	400 lbs.	24"	19"	18"	Enhanced 3-Position	No	No / No	USB-A & USB-C	Cup holder; arm storage; magnetic hand-control docking; integrated battery backup; fully padded chaise	17"
Golden MaxiComforter PR535	Petite/Small	Not specified	375 lbs.	21"	19"	20"	MaxiComfort 2 Comfort Zones	Optional Heat & Vibrating Massage	No / No	No	AutoDrive hand control; signature seat; pocketed coil seat; SmartTek diagnostic system; stain defense; optional footrest extension; 9V battery backup	22"
Pride VivaLift! Radiance	Petite Wide	5'4" & below	400 lbs.	24"	18"	19.5"	True-Infinite 4-Motor	3-Zone Heat No massage listed	Yes / Yes	USB + wireless charging	Power backrest; standard footrest extension; 2 memory positions; lithium battery backup; fully padded chaise	14"

Note: Specifications shown are for the petite or petite-wide model of each lift chair. Golden's sell sheets do not specify a recommended user height for the PR531 Petite/Small or PR535 Petite/Small. Measurements and features may vary by model and configuration.

FOR FREEDOM MOBILITY SOLUTIONS EXCLUSIVE USE ONLY