

Space-Saving Lift Chair Comparison Chart

Freedom Mobility Solutions | Compact Lift Chair Models for Smaller Spaces

Model	Model #	Sizes	Overall Width	Wall Clearance	Weight Capacity	Positioning	Heat / Massage	Power Headrest / Lumbar	Other Features
Pride Genesis	LC-150	Small / Medium	32-33 in.	17-18.5 in.	325 lbs.	Traditional 3-Position	No	No / No	Compact overall width; fully padded chaise; battery backup
Pride VivaLift! Metro 2	PLR-925M	Medium	33 in.	14 in.	400 lbs.	True-Infinite	No	Yes / Yes	Power headrest/lumbar; footrest extension; 2 memory positions; USB
Golden EZ Sleeper Slim	PR750	Medium	33.5 in.	21.5 in.	400 lbs.	4-Zone MaxiComfort	Optional	Yes / Yes	Slim track-arm design; USB; lithium backup; ZG+
Golden EZ Sleeper Slim w/ Twilight	PR764	Small / Medium	33.5 in.	21.5 in.	400 lbs.	5-Zone MaxiComfort + Twilight	Optional	Yes / Yes	Twilight Tilt; USB; lithium backup; ZG+
Golden Comforter	PR531	7 sizes	Varies by size*	Varies by size*	400-500 lbs.	3-Position / 1 Zone	Optional	No / No	Petite through wide/tall sizing; USB hand control; optional footrest extension

Note: Space-saving does not mean every chair is a wall-hugger. These models were selected for compact overall width, reduced wall clearance, or configurations that work well where floor space is limited. *Golden Comforter dimensions vary by size. For example, Golden currently lists a 34 in. overall width and 16 in. wall clearance for the displayed configuration. Always verify the selected size before purchase.

FOR FREEDOM MOBILITY SOLUTIONS EXCLUSIVE USE ONLY