



3 Motor Full-Electric Bed

WITH INCLUDED LOW BED CAPABILITY

The **ProBasics Full Electric Bed With Included Low Bed Capability** combines the features of a standard full-electric bed with the safety and convenience of a low bed, where rails are not desirable or appropriate, but potential injuries from falling out of bed are still a concern. Dual-purpose bed ends with telescoping legs accommodate multiple bed height profiles from low to standard to high and are self-leveling.

The split-spring design incorporated into the low bed profile still allows for the easy, one-person delivery but with all the benefits of a low bed.



Lowest Height Position for Low Bed at 10.3" -23" Standard Height Position 12.7" min - 25.3" max

FEATURES

- Bed ends convert from low bed to conventional full electric with telescoping legs and are self-leveling
- Strong with a welded steel split spring design and reinforced head section
- Easy assembly with motors & electronics fully removable for easy cleaning of frame
- The easy-to-use pendant control provides motorized positioning of the upper body and knees, and adjusts the bed frame height.
- The head and foot sections can be raised or lowered both simultaneously or independently
- Composite head and footboard panels are both strong and durable
- In case of a power failure head and foot section can be lowered by adding two 9V batteries on to the junction box via a simple connection
- HCPCS Code: E0297*

PRODUCT SPECIFICATIONS

Bed Height:	10.3" min. – 25.3" max.
Sleep Surface:	80" L x 36" W
Overall Bed:	88" L x 36" W
Weight Capacity:	450 lb
Maximum Patient Weight:	350 lb
Limited Warranty:	5 year on frame 2 year on motor - mechanical/electrical

SHIPPING SPECIFICATIONS

Component:	Carton Size	UPC
3MHBFESSLO	45.5" X 10.5" X 37"	815067078615
3MHBFEVELO	38" X 7.5" X 32.25"	815067078622
Shipping Weight:	185.4 lbs	

Item #	Description	UOM
PBFEB3MLO	ProBasics 3 Motor Full Electric Bed with Included Low Bed Capability	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice. Information is accurate at the time of printing, but specifications are subject to change.