

## Viper Plus GT Wheelchair

HCPCS: K0004-1,K0004+K0195+E0990



FIGURE A



FIGURE B



FIGURE C

- Built-in seat rail extensions and extendable upholstery easily adjusts seat depth from 16" to 18" (Patent pending) (Figure A - C)
- All-Aluminum frame
- Dual Axle allows for multiple seat-to-floor height positions
- Adjustable height, flip-back, removable, padded arms
- Weighs under 35 lbs in basic configuration (excluding front rigging)
- Over-center cross brace eliminates need for seat guides Works with seating and positioning accessories
- Adjustable angle back easily adjusts from 5° - 20°
- Adjustable height back adjusts from 17" - 19" in 1" increments, adaptable for individual patient needs
- Durable, heavy-gauge flame-retardant nylon upholstery with hook and loop fastener flaps
- Quick release, 24" composite Mag-style wheels are lightweight and maintenance free
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Adjustable angle caster forks have three height adjustments and angle adjustability
- Tool-free adjustable swing-away footrests also feature easy-to-access release handles
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight
- Adjustable anti-tipper with flip-up wheels Never needs to be removed
- Comes with push-to-lock wheel locks
- Patent Pending

### PRODUCTS

#### Description

PLA418FBDAARAD-ELR	ViperPlsGT,18"FlpBkAdjHtDskArm,EL,R1/ea
--------------------	---

### SPECIFICATIONS

#### Description

Warranty	Limited Lifetime
Weight Capacity	300 lbs   350 lbs
Height	36"
Seat To Floor Height	17.5" - 19.5"
Carton Shipping Weight	51 lbs   52 lbs   53 lbs   54 lbs
Arm Rest Type (s)	Flip Back, Detachable, Adj. Height Desk Arm   Flip Back, Detachable, Adj. Height Full Arm   Flip-Back, Detachable, Adj. Height Desk Arm
Back Height	17" - 19"

### ACCESSORIES

#### Description

STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA

ACCESSORIES	Description
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS807	Anti-Tips,Adj.w/Whls,Blk,1pr/bx,10/cv, BX
STDS828GT	AmputeeAdaptersViperGT,1pr, PR
STDS8HK	Super Hemi Kit Viper / Viper GT, 1kit/cs, CS
STDS4056	One Arm drive for ViperGT, EA
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
PARTS	Description
STDS2369BLK	Back Uphol 22"x18",Black,1/bx, BX
STDS1A3510	Rear Wheel,VIPER GT,1/EA, EA
STDS1A3502	Frame Assist, 22" ViperPlusGT, 1/ea, EA
STDS1A3511	Quick Release Rear Axle, Viper GT, 1/ea, EA
STDS5D1012	24" Rear Wheel Bearing, C4,1 PR, PR
STDS1073	Spacers & Screws for Handrim, 8/set, EA
STDS2359BLK	Back Uphol 18"x16",Black,1/bx, BX
STDS2361BLK	Back Uphol 20"x18",Black,1/bx, BX
STDS1A3557	Backrest Cane,VIPER GT,1PR, PR
STDS2A3250	BackAngle/Height Adj.Vip.Plus,1/ea, EA
STDS2363BLK	Back Uphol 16"x18",Black,1/bx, BX
STDS1A3513	Front Axle Bolt , 1/pr, PR
STDS2A3511	Hand Rim,Viper GT,1ea, EA
STDS2364BLK	Seat Uphol 18"x16",Black,1/bx, BX
STDS2366BLK	Seat Uphol 20"x18",Black,1/bx, BX
STDS2368BLK	Seat Uphol 22"x18",ViperGT,1/ea, BX
STDS2356BLK-AD	Seat UpholsteryNylon,16"ViperPlusGT,1ea, EA
STDS2364BLK-AD	Seat UpholsteryNylon,18"ViperPlusGT,1ea, EA
STDS1A35DL	Adj.DeskArm,Lt.Vip.GT,1/ea, EA
STDS1A35DR	Adj.DeskArm,Rt.Vip.GT,1/ea, EA